
SteriMax

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Instruction manual EN

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Impfösen - Sterilisator

Inoculation Loop - Sterilizer

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
Fig. C




Fig. D

- 1 - Master switch
- 2 - Residual heat display
- 3 - Annealing tube consists of:
 - 3a: outer tube: Quartz glass heating element
 - 3b: inner tube: Ceramic annealing tube
- 4 - Handle for angular adjustments
- 5 - Power socket
- 6 - Fuse holder
- 7 - Splash guard
- 8 - Depth stop made of quartz glass (inside the ceramic annealing tube 3b)
- 9 - Retaining screw for splash guard
- 10 - Connector
- 11 - Power cable
- 12 - Power plug (plug type may vary for different countries)
- 13 - Tray

Illustration of the 7 angular adjustments

Read these operating instructions carefully to familiarize yourself with the product before using. Keep the manual in a safe place for future reference.

Table of contents

The range.....18

Use.....18

General note.....18

General safety advise.....18

Warnings.....19

1. Initial operation.....21

2. Turning the device On / Off.....21

3. Residual heat display.....21

4. Sterilising instruments.....21

5. Angle adjustment.....21

6. Cleaning and replacing the quartz glass / ceramic annealing tube.....22

7. Replacing the fuse.....22

8. Holder installation (optional).....23

9. Warranty.....23

10. Troubleshooting guide.....23

11. Maintenance instructions.....24

 11.1 Replacing the quartz glass heating element.....24

 11.2 Removal of the the quartz glass heating element.....24

 11.3 Installation of the the quartz glass heating element.....25

Technical Data (5.001.000, 5.002.000, 5.003.000).....26

Technical Data (5.101.000, 5.102.000, 5.103.000).....27

Accessories.....28

Spare Parts.....28

Declaration of Conformity.....29

The range

SteriMax - Inoculation Loop Sterilizer:

Length 118 mm, Ø 15 mm:

230 V / AC Art. No.: 5.001.000

120 V / AC Art. No.: 5.002.000

100 V / AC Art. No.: 5.003.000


Length 161 mm, Ø 20 mm:

230 V / AC Art. Nr.: 5.101.000

120 V / AC Art. Nr.: 5.102.000

100 V / AC Art. Nr.: 5.103.000

- with 7 angular adjustments
- 1 screwdriver (R1)
- 1 inoculation loop holder (R2)
- 1 inoculation loop Ø 3 mm (R3)
- Tray for 3 WLD-TEC loop holder (13)
- Power cable (11)
- 2-year warranty

Use

The inoculation loop sterilizer Sterimax is ideal for laboratory use for sterilizing inoculation loops, tweezers and micro-instruments made of heat resistant materials such as stainless steel or platinum-iridium alloys.

General note


Ceramic insulated heating element: The initial brown color of the heating element insulation is normal. The heating element decolors during use.

General safety advise

Only operate the annealing system under constant supervision!

Clean the ceramic annealing tube frequently and check for damage. If the ceramic annealing tube (3b) exhibits any cracks, chips or similar, do not use the device! Replace the ceramic annealing tube. Check also if the quartz glass heating element (3a) is damaged! For further information see paragraph 6 of this manual.

- All users who have been assigned to use this device must read and understand this manual or have been instructed by a competent person in such a way that they can use this device without causing danger.
- When unpacking, check for possible shipping damage and do not operate if there are visible damages.
- After using the device or if it is unused for longer periods, turn off the power at the master switch (1).
- Do not operate the device near flammable liquids, materials or in potentially explosive areas.
- When using the loop holder, avoid leaving the loops and needles within the annealing tube (3) for extended periods of time. Doing so will overheat the loop rest causing a potential safety risk and may damage your loop or needle. Furthermore melted parts can damage the annealing tube and are a safety risk.
- Never insert disposable plastic loops or needles into the annealing tube (3) with a temperature resistance less than 1200°C.
- Never use the device without insert the ceramic annealing tube (3b).
- Using a two-pin AC adapter / connecting cable is NOT recommended.
- Operatin temperature 5°C bis 40°C. Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.

Warning note in order to avoid the risk of electric shock and / or possible death.

- Make sure that the unit does not come in contact with water or other liquids as long as the plug is in the power socket, especially if you use it in the laboratory near water or liquids.

- Make sure that during operation the power cable (10, 11, 12) never gets wet or damp. Locate the cable so that it is not trapped, jammed or otherwise damaged. If the power cable (10, 11, 12) is damaged, the cable must first be replaced before using this device.

- After using the device always unplug the power plug (12) from power outlet to avoid unintentional switching on.

- During malfunctions or when cleaning the device, always disconnect the power plug (12).

- If the quartz glass heating element (3a) exhibits any cracks, chips or similar, turn off the device and unplug the mains plug (12) from power socket. Do not use the device! Replace the heating element before using the device again.
- If the ceramic annealing tube (3b) exhibits any cracks, chips or similar, turn off the device and unplug the mains plug (12) from power supply (5). Do not use the device! Replace the ceramic annealing tube (see paragraph 6). Check also if the quartz glass heating element (3a) is damaged.

Warning note in order to avoid the risk of fire and/or injuries.

- During the operation of the device, its components may get extremely hot. Therefore, only touch the controls or the handle (4) for the adjustment of the device's angular settings.

- Never place your hands or any other body part above the opening of the annealing tube (3)!
- Combustible materials may ignite within the immediate area of or inside the annealing tube (3)! Never place the device close to combustible objects and / or materials. Never insert combustible materials or liquids into the annealing tube (3).
- Never insert instruments or objects where flammable liquids or material adhere to these instruments or objects into the annealing tube (3).
- Never cover up the annealing system as long as it is in use or its casing is still hot.
- Operate the unit only on a stable, non-slip and level surface.

- During the heating-up of the annealing tube the display of the residual heat (2) will be delayed. It is possible that the tube's casing is already hot during the heating phase or before the residual heat display is turned on.


1. Initial operation

Please handle this sterilizer always with care as described in the following points of this manual.


It is essential to operate the annealing system only on an even and level surface. Never operate the device next to highly inflammable liquids and / or materials, or inside enclosed rooms which can be subject to explosion hazards.

Connect the power cable (11) with the connector (10) to the device's power socket (5). Insert the power plug (12) into a live power outlet suitable for your device (230V AC, 120V AC, 100V AC).

2. Turning the device On / Off

The device may be turned off and on, using the master switch (1). As soon as the device is turned on, a green light will illuminate the master switch (1). After the device has been turned on, the annealing tube (3) will require several minutes to heat up to a temperature of 900 - 950°C.

3. Residual heat display

The red residual heat display (2) warns that the casing of the annealing tube is hot. The residual heat display (2) remains illuminated after the device has been switched off and until the casing of the annealing tube has been cooled off.

During the heating-up of the annealing tube, the display of the residual heat (2) will be delayed. It is possible that the tube's housing is already hot during the heating phase or before the residual heat display (2) is turned on. If the power supply to the device is interrupted, the residual heat display (2) will extinguish, even if the casing of the tube is still hot!

4. Sterilising instruments

Insert the instrument carefully (e.g. loop with holder) in the annealing tube (3). Put it at least half way inside the annealing tube to avoid evaporation of contaminated residues. The instrument will be sterilized in about 5 - 7 seconds* and it will start to glow. Not in use instruments can be placed on the tray (13).

During the operation of the device, its components may get extremely hot. Therefore, only touch the controls or the handle (4) for the adjustment of the device's angular settings. Never place your hands or any other body part above the opening of the annealing tube (3)!

5. Angle adjustment

The device is equipped with an angle adjustment mechanism that can lock itself into 7 different positions. The angles can be easily set, simply adjust with the handle (4) on the casing of the annealing tube. The various angle settings are shown in the figure on page 16.

* tested with inoculation loops WLD-TEC (Art.-No.: 6.000.373)

6. Cleaning and replacing the quartz glass / ceramic annealing tube

Before cleaning / replacing, wait until the device cooled down, then remove the power plug (12) from the power outlet and remove the connector (10) from the power socket (5).

For cleaning / replacing purposes the ceramic annealing tube (3b) can be removed easily at the rear side of the unit. For this purpose unscrew the retaining screws (9) of the splash guard (7) on the back of the device. Remove splash guard (7) and take out depth stop (8). Now pull out carefully the ceramic annealing tube (3b). The annealing tube is made of ceramic and can be easily cleaned while slightly wet, using a commercially available non-alkali cleaning agent. While the ceramic annealing tube (3b) is removed, the heating element made of quartz glass (3a) can be easily cleaned as well, slightly wet, using a commercially available non-alkali cleaning agent. Pipe cleaners or bottle brushes are recommended for the cleaning. After the cleaning, put the ceramic annealing tube (3b) and the depth stop (8) back into the quartz glass heating element (3a) and reinstall the splash guard (7).

Do not heat up the unit without assembled ceramic annealing tube (3b)!


Fig. 6/1


Fig. 6/2


Fig. 6/3

7. Replacing the fuse

When replacing the fuse, first remove the power plug (12) from the power outlet and remove the connector (10) from power socket (5).

The device is equipped with a fuse (Flink). Always replace a defective fuse! To do this, pull out the fuse holder (6). (For easy removal of the fuse holder, use a suitable tool and carefully pull on the small tab, intended for this purpose). Replace the blown fuse with a new one and push the fuse holder (6) back into its original position. A spare fuse has been provided and can be found in the retainer (6a) of the fuse holder (6).


Fig. 7/1

Fuses

Sterimax 230 V (5.001.000): Fuse 800 mA (F)
Sterimax 120 V (5.002.000): Fuse 1,6 A (F)
Sterimax 100 V (5.003.000): Fuse 2,0 A (F)
Sterimax 230 V (5.101.000): Fuse 1,25 A (F)
Sterimax 120 V (5.102.000): Fuse 2,0 A (F)
Sterimax 100 V (5.103.000): Fuse 2,5 A (F)


Fig. 7/2


8. Holder installation (optional)

Before installing the holder turn off the device, wait for it to cool down and unplug it!

Two optional holders are available for the SteriMax:

1. Holding device, specially developed for WLD-TEC inoculation loop holders, made of stainless steel, available for all angle settings, **Art. No.: 5.000.100.**
2. Support for all inoculation loop holders, made of stainless steel; can only be used in the first three angle settings, **Art. No.: 5.000.101.**

EN

Both holders will be attached to the underside of the annealing tube's casing using two short screws (included with delivery of the holder). When installing the holders, first push the screws (**Z1**) through the slots of the holder (**Z2**) and slightly tighten the screws in the tapped hole on the underside of the tube's casing. Subsequently, insert an inoculation loop holder and slide the holder into the desired position. Tighten the screws (**Z1**).

9. Warranty

All SteriMax are covered under our two-year manufacturer warranty against any manufacture defects in material and workmanship, except parts to wear and tear (Quartz glass heating element **3a** & ceramic annealing tube **3b**).

The WLD-TEC warranty guarantees all SteriMax under normal usage conditions and does not cover any damages as a direct result of user misuse or / and abuse.

The warranty is void upon any unauthorized servicing, disassembly or modifications.

10. Troubleshooting guide

• Device does not work

Ensure the mains cable is connected properly.
(see paragraph 1).

Check fuse, replace if necessary.
(see paragraph 7).

• Annealing tube (3) does not heat up, even though the device is turned on and the master switch's green light is illuminated.

Heating element may be defect. Have it replaced by qualified person only! (see paragraph 11).

11. Maintenance instructions

This section describes how to replace the quartz glass heating element.

The removal and installation of the quartz glass heating element may only be carried out by qualified personnel!

11.1 Replacing the quartz glass heating element

The long-lasting quartz glass heating element is to a large extent maintenance and wear-free. However, when working over a longer period at higher voltages, abrasive substances may dirty or damage the white insulation material of the heating element; in this case, replacement may be required.

For safety reasons the entire quartz glass heating element has to be replaced if the quartz glass heating element (3a) becomes damaged and cracks, chips or similar damages are noticeable.


Fig. 11/1


Fig. 11/2

11.2 Removal of the the quartz glass heating element

Before removing the heating element, wait until the device cooled down, then remove the power plug (12) from the power outlet and remove the connector (10) from the power socket (5). Danger of electric shock!

Remove retaining screws (9) of the splash guard (7) on the back of the device. Remove splash guard (7) and pull out the depth stop (8) and the ceramic annealing tube (3b) (see paragraph 6). Subsequently, remove both screws (E1) from the casing and remove carefully (E2) the rear panel.


Fig. 11/3

Pull the heating element (3a) just a few centimetres out of the casing. While using a suitable screwdriver, loosen both screws (E4) on the power side of the ceramic terminal (E5). Pull both black wires (E6) from the ceramic terminal (E5). Completely remove the heating element (3a) from the casing.


Fig. 11/4

11.3 Installation of the the quartz glass heating element

Carefully insert the new heating element **(3a)** into the casing. (Heating element is fragile) Insert both black wires **(E6)** into the ceramic terminal **(E5)** of the new heating element **(3a)**. Use a suitable screwdriver to firmly tighten both screws **(E4)** of the new heating element. The black wiring shall be routed such as to maintain a minimum distance of 1 - 2 centimetres from the white insulation of the heating element.


Fig. 11/5

When tightening the screws, ensure to hold the connector piece (E5) in order to prevent it from twisting. Risk of short circuiting!

Now, push the heating element with front end of the quartz glass tube into the front guide **(E7)** of the casing. Re-install the rear panel of the casing **(E2)**. Ensure to slide the rear guide **(E8)** over the rear end of the quartz glass tube. At the same time, adjust the casing's wall **(E2)** to allow for the proper attachment of the casing's screws **(E1)**. Tighten the casing's screws **(E1)**.

Do not pinch the wiring between the casing's wall (E2) and/or other components of the casing! Do not pull on the connection wire (E3) of the residual heat sensor and do not damage the insulation. Danger of electric shock! If the wires of the residual (E3) have moved, place them again close to the lower casing part.

Finish the installation by re-inserting the ceramic annealing tube **(3b)** and depth stop **(8)** into the quartz glass heating element **(3a)** and attach the splash guard **(7)** using both screws **(9)**.

Do not heat up the unit without assembled ceramic annealing tube (3b)!

Service address:

WLD-TEC GmbH
Production & Service
Halle-Kasseler-Str.49
D-37318 Arenshausen
Germany

Phone: +49 36081 68940
Fax: +49 36081 68942
Email: sales@wld-tec.com
Internet: www.wld-tec.com



Technical Data (5.001.000, 5.002.000, 5.003.000)

Mecanical

Quartz glass heating element:	Ø 21 mm
Ceramic annealing tube:	Ø 15 mm
Ceramic annealing tube (length):	118 mm
Angular adjustments:	7
Weight:	960 g
Measurements (w x h x d):	175 x 129 x 139 mm
Casing:	stainless steel
Heating element:	quartz glass
Annealing tube:	ceramic, 60% Al ₂ O ₃
Heating element Insulation:	MgO 19-26%, SiO ₂ 72-77%

Electical

Sterilization temperature:	900°C - 950°C (1650°F - 1742°F)
Power supply:	Art. No.: 5.001.000: 220 - 240 V / 50/60 Hz / 160 W Art. No.: 5.002.000: 112 - 128 V / 50/60 Hz / 160 W Art. No.: 5.003.000: 95 - 105 V / 50/60 Hz / 160 W
Fuses:	Art. No.: 5.001.000: 800 mA (F) Art. No.: 5.002.000: 1,6 A (F) Art. No.: 5.003.000: 2,0 A (F)
Protection class:	I
Heating up time:	700°C / 11 minutes; 800°C / 14 minutes
Sterilization time:	5 - 7 seconds*
Heat transfer:	a stainless steel surface, located above, warms up at motionless air and at an ambient temperature of 20°C; Distance: 20 cm: 65°C, 30 cm: 34°C, 40 cm: 28°C

Enviromental conditions

operating temperature 5°C bis 40°C. Maximum relative humidity 80 % for temperatures up to 31°C, decreasing linearly to 50 % relative humidity at 40°C

Safety features

Residual heat display: indicates a hot annealing tube

Licences

CE: EN 61326-1, EN 61010-1, EN 61010-2-010
EU guidelines: 2014/30/EU, 2014/35/EU, 2011/65/EU

*tested with inoculation loops, WLD-TEC (Art.-No.: 6.000.373)

Technical Data (5.101.000, 5.102.000, 5.103.000)

Mecanical


Quartz glass heating element:	Ø 26,6 mm
Ceramic annealing tube:	Ø 20 mm
Ceramic annealing tube (length):	161 mm
Angular adjustments:	7
Weight:	1060 g
Measurements (w x h x d):	175 x 129 x 169 mm
Casing:	stainless steel
Heating element:	quartz glass
Annealing tube:	ceramic, 60% Al ₂ O ₃
Heating element Insulation:	MgO 19-26%, SiO ₂ 72-77%

EN

Electical

Sterilization temperature:	900°C - 950°C (1650°F - 1742°F)
Power supply:	Art. No.: 5.101.000: 220 - 240 V / 50/60 Hz / 225 W Art. No.: 5.102.000: 112 - 128 V / 50/60 Hz / 225 W Art. No.: 5.103.000: 95 - 105 V / 50/60 Hz / 225 W
Fuses:	Art. No.: 5.101.000: 1,25 A (F) Art. No.: 5.102.000: 2,0 A (F) Art. No.: 5.103.000: 2,5 A (F)

Protection class:	I
Heating up time:	700°C / 15 minutes; 800°C / 22 minutes
Sterilization time:	5 - 7 seconds*
Heat transfer:	a stainless steel surface, located above, warms up at motionless air and at an ambient temperature of 20°C; Distance: 20 cm: 65°C, 30 cm: 40°C, 40 cm: 32°C

Enviromental conditions

operating temperature 5°C bis 40°C. Maximum relative humidity 80 % for temperatures up to 31°C, decreasing linearly to 50 % relative humidity at 40°C

Safety features

Residual heat display:	indicates a hot annealing tube
------------------------	--------------------------------

Licences

CE:	EN 61326-1, EN 61010-1, EN 61010-2-010
EU guidelines:	2014/30/EU, 2014/35/EU, 2011/65/EU

*tested with inoculation loops, WLD-TEC (Art.-No.: 6.000.373)

Accessories

Holding device (for WLD-TEC loop holder):	Art. No.: 5.000.100
Universal support (for all loop holder):	Art. No.: 5.000.101
Tray (for left side, stainless steel):	Art. No.: 6.000.341
Loop holder (stainless steel):	
length 245 mm	Art. No.: 6.000.360
length 215 mm	Art. No.: 6.000.365
Inoculation Loops (special stainless steel):	
Ø 3 mm:	Art. No.: 6.000.373
Ø 5 mm:	Art. No.: 6.000.375
Inoculation Loops (Platin / Iridium):	
Ø 1 mm:	Art. No.: 6.000.381
Ø 3 mm:	Art. No.: 6.000.383
Ø 5 mm:	Art. No.: 6.000.385

Spare Parts

Quartz glass heating element 230 V (for 5.001.000):	Art. No.: 5.000.202
Quartz glass heating element 120 V (for 5.002.000):	Art. No.: 5.000.203
Quartz glass heating element 100 V (for 5.003.000):	Art. No.: 5.000.204
Quartz glass heating element 230 V (for 5.101.000):	Art. No.: 5.100.202
Quartz glass heating element 120 V (for 5.102.000):	Art. No.: 5.100.203
Quartz glass heating element 100 V (for 5.103.000):	Art. No.: 5.100.204
Ceramic annealing tube (for 5.001.000, 5.002.000, 5.003.000):	Art. No.: 5.000.205
Ceramic annealing tube (for 5.101.000, 5.102.000, 5.103.000):	Art. No.: 5.100.205
Fuse 230 V (800 mA, F) (for 5.001.000):	Art. No.: 5.001.300
Fuse 120 V (1,6 A, F) (for 5.002.000):	Art. No.: 5.002.300
Fuse 100 V (2,0 A, F) (for 5.003.000):	Art. No.: 5.003.300
Fuse 230 V (1,25 A, F) (for 5.101.000):	Art. No.: 5.101.300
Fuse 120 V (2,0 A, F) (for 5.102.000):	Art. No.: 5.102.300
Fuse 100 V (2,5 A, F) (for 5.103.000):	Art. No.: 5.103.300

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