

TECHNICAL DATA SHEET

Article no. : 8485.0500

MAXIMUM RECOVERY DILUENT (MRD)

ALSO KNOWN AS

UNIVERSAL DILUENT

SPECIFICATION

Isotonic diluent for the maximal recovery of stressed microorganisms according to ISO standards.

FORMULA * IN G/L

pH final 7,0 ±0,2 at 25 °C *Adjusted and /or supplemented as required to meet performance criteria

DIRECTIONS

Dissolve 9.5 g of powder in 1 L of distilled water and distribute into suitable containers. Sterilize in the autoclave at 121°C for 15 minutes.

DESCRIPTION

This formulation combines the osmotic pressure of the physiological saline solution with the protective action of the peptone to obtain good recovery of stressed microorganisms.

The sodium chloride ensures isotonic conditions and the low concentration of peptone does not allow cellular growth in the short period (2-3 hours) of time required for the preparation of the dilution bank of the sample.

TECHNIQUE

According to the ISO method, the sample is diluted in a ratio 1:10 with the Maximum Recovery Diluent and homogenized by a vortex mixer or Stomacher®. After a short period (10-15 minutes) of rest, a 1/10 dilution bank with the same diluent is prepared following standard procedures. Plates are inoculated using the range of different concentrations.

QUALITY CONTROL

 Incubation Temperature:
 35°C ±2,0

 Incubation Time:
 Recovery 24 h

 Inoculum:
 1000-10000 CFU / tube (Productivity) at T0, 45 minutes and 1 h. (20-25°C) ; according ISO 11133:2014/Amd 1:2018

 Microorganism
 Growth
 Remarks

-		
Staphylococcus aureus ATCC® 6538	Good	Recovery ±30% T0, in TSA
Pseudomonas aeruginosa ATCC® 9027	Good	Recovery ±30% T0, in TSA
Escherichia coli ATCC [®] 8739	Good	Recovery ±30% T0, in TSA
Candida albicans ATCC [®] 10231	Good	Recovery ±30% T0, in SDA



Th. Geyer GmbH & Co. KG Domierstr. 4 – 6 D-71272 Renningen Tel: +44 7159 1837-0 Fax: +49 7159 1837-710 renningen@thgeyer.de www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600) IBAN DE856005101000203802 Postbank Stuttgart (Swift/BIC PBNKDEFFXXX) IBAN DE32600100700000020708 Deutsche Bank (Swift/BIC DEUTDESSXXX) IBAN DE06600700700125518100 St.-Nr. 70093/40018 / USt-IdNr. DE147510304 Amtsgericht Stuttgarf / HRA-Nr. 254140 Persönlich haftende Gesellschaftrein: Geyer Beteiligungsgesellschaft mbH Amtsgericht Stuttgarf / HRB-Nr. 252035 Geschäftschner. Lutz-Alexander Geyer / Oliver-Alexander Geyer / André Meise / Ralf Streicher



REFERENCES

- · ISO 6887-1: 1999 Microbiology of food and animal feeding stuffs. Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: General rules for the preparation of the initial suspension and decimal dilutions - Part 2 (2003): Specific rules for the preparation of meat and meat products
- · ISO Standard 8199 (2018) Water Quality General requirements and guidance for microbiological examinations by culture.
- · ISO 8261: 2001 Standard. Milk and milk products General guidance for the preparation of test samples, initial suspension and decimal dilution for microbiological examination.
- . ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- · ISO 16212 Standard (2017) Cosmetics Microbiology Enumeration of yeast and mould.
- · ISO 18415 Standard (2017) Cosmetics Microbiology Detection of specified and non-specified microorganisms.
- · ISO 18416 Standard (2015) Cosmetics Microbiology Detection of Candida albicans.
 · ISO 21149 Standard (2017) Cosmetics Microbiology Enumeration and detection of aerobic mesophilic bacteria.
- · ISO 21150 Standard (2015) Cosmetics Microbiology Detection of Escherichia coli.
- · ISO 22717 Standard (2015) Cosmetics Microbiology Detection of Pseudomonas aeuruginosa.
- · ISO 22718 Standard (2015). Cosmetics Microbiology Detection of Staphylococcus aureus.

STORAGE

For laboratory use only. Keep tightly closed, away from bright light, in a cool dry place (+4 °C to 30 °C).



Th. Geyer GmbH & Co. KG Dornierstr. 4 - 6 D-71272 Renningen Tel.: +49 7159 1637-0 Fax: +49 7159 1637-7 -710 renningen@thgeyer.de www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600) IBAN DE85600501010002036302 Postbank Stuttgart (Swift/BIC PENKDEFFXXX) IBAN DE32600100700000020708 Deutsche Bank (Swift/BIC DEUTDESSXXX) IBAN DE06600700700125518100 St.-Nr. 70093/40018 / USt-IdNr. DE147510304 Amtsgericht Stuttgart / HRA-Nr. 254140 Persönlich haftende Gesellschafterin: Gever Betelligungsgesellischaft mbH Amtsgericht Stuttgart / HRB-Nr. 252035 Geschäftsführer: Lutz-Alexander Gever / Oliver-Alexander Gever / André Meise / Ralf Streicher