

Technical Datasheet

Article No. 1196

Palladium matrix modifier

10000 µg/ml in 15% HNO₃

For laboratory use.

Parameter	Value
Appearance	Clear, yellowish liquid

All standard components have been pre-qualified/verified before use. All analytical instruments and instrumentation were pre-calibrated. The actual concentrations given below are based on this preparation method and the impurities of the compounds.

Raw material	Palladium
Purity of Raw Material	min. 99.99 %
Assay mg/kg	min. 10.000
Assay mg/l	min. 9.8000