

TECHNICAL DATA SHEET

Article no. : 8813.0500

LISTERIA SELECTIVE AGAR BASE ACC. TO OTTAVIANI AND AGOSTI

SPECIFICATION

Medium for the isolation of *Listeria spp.* and the presumptive identification of *L. monocytogenes* acc. to ISO 11290-1 and 11290-2 standard.

FORMULA * IN G/L

Meat peptone.....	18.000	Lithium chloride.....	10.000
Tryptone.....	6.000	Disodium phosphate anhydrous.....	2.500
Yeast extract.....	10.000	5-bromo-4-chloro-3-indolyl-	
Sodium pyruvate.....	2.000	β-D-glucopyranoside.....	0.050
Dextrose.....	2.000	Agar.....	12.000
Magnesium glycerophosphate.....	1.000		
Magnesium sulphate.....	0.500	Final pH 7.2 ± 0.2 at 25 °C	
Sodium chloride.....	5.000		

*Adjusted and /or supplemented as required to meet performance criteria

DIRECTIONS

Suspend 33 g of powder in 476 ml of purified water and bring to the boil with constant stirring. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 45-50 °C and add 1 bottle of Listeria enrichment supplement Ottaviani & Agosti (Art. 8815.0010) and 1 vial of Listeria selective supplement Ottaviani & Agosti (Ref. Art. 8814.0010). Homogenize by mixing and distribute in Petri dishes. The solidified medium appears turbid.



DESCRIPTION

Completed with all its supplements the Agar Listeria O&A is a selective and differential medium for the detection of *Listeria species* and the presumptive identification of *Listeria monocytogenes*.

The selectivity is achieved by the high concentration of lithium chloride and the mixture of antimicrobics. The differential activity is due to the chromogenic substrate to detect the β -glucosidase enzyme that is present in all *Listeria species*.

The specific identification is obtained by the L- α -phosphatidylinositol, that acts as substrate for a phospholipase C that is present only in *Listeria monocytogenes* and some strains of *Listeria ivanovii*.

The combination of both substrates allows the differentiation *L. monocytogenes* that produces colonies blue-green in colour but surrounded by an opaque zone from the other *Listeria species* that growth with blue-green colonies without any halo. This differentiation is evident after incubate the plates for 24 \pm 2 hours at 37 °C.

Sometimes, especially with highly contaminated samples it is possible that can growth some colonies, white in colour, that are not Listeria. In this case an enrichment step is recommended prior to previous to the plate inoculation.

Observations: Most *Listeria ivanovii* also produce an opaque halo around the colonies after 48 h of incubation. This presumptive evidence must be confirmed by performing the biochemical or serological identification tests (Rhamnose / Xylose sugar fermentation, hemolysis tests, CAMP test, etc.) or any test confirming the species without hesitation.

Remarks:

Enrichment Supplement for Agar Listeria O&A (335032ZF):
1 vial sufficient amount for 500 ml complete medium

L- α -phosphatidylinositol..... 1.00 g
Sterile Distilled Water..... 24.00 ml

Selective Supplement for Agar Listeria O&A (8814.0010):
1 vial sufficient amount for 500 ml complete medium

Nalidixic acid..... 10 mg
Ceftazidime..... 10 mg
Cycloheximide..... 25 mg
Polymyxin B sulphate..... 38350 IU

Reconstitute the original freeze-dried vial by adding 1 vial with
Sterile distilled water..... 6.00 ml

TECHNIQUE

There are many standardised methodologies (ISO, FDA-BAM, AOAC, AFNOR, etc.). The technician must follow the protocol validated in his laboratory.

QUALITY CONTROL

Incubation Temperature: 37 \pm 1.0 °C

Incubation Time: 44 \pm 4 h

Inoculum: Practical range 100 \pm 20 CFU. min. 50 CFU (productivity)/ 10⁴-10⁶ CFU (selectivity)/ 10³-10⁴ CFU (specificity), according to ISO 11133:2014/Amd 1:2018.

Microorganism	Growth	Remarks
<i>Escherichia coli</i> ATCC® 25922	Inhibited	-
<i>Listeria monocytogenes</i> ATCC® 13932	Productivity > 0.50	green-blue colonies surrounded by opaque halo
<i>Listeria monocytogenes</i> ATCC® 35152	Productivity > 0.50	green-blue colonies surrounded by opaque halo
<i>Listeria innocua</i> ATCC® 33090	Good (specificity)	green-blue colonies without opaque halo
<i>Enterococcus faecalis</i> ATCC® 29212	Inhibited	-



REFERENCES

- Artault, S., J.L. Bind, Y. Delaval, N. Dureuil, N. Gallart (2000) AFNOR validation of the ALOA method for the detection of *Listeria monocytogenes* in foodstuffs. Coll. Soc. Fran. Microbiol. 19-20 Oct. Paris.
- Bannerman, E.S. & J. Bille (1988) A new selective medium for isolating *Listeria* from heavily contaminated material. Appl. Environ. Microbiol. 54:1:165-167.
- Greenwood, M., C. Willis, P. Dosweell, G. Allen & K. Pathak (2005) Evaluation of chromogenic media for the detection of *Listeria* species in food.
- Hitchins, A.D. & K. Jinneman (1998) *Listeria monocytogenes* in FDA-BAM 8th edition Revision A. Updater January 2003. AOAC Intl. Gathersburg. MD. USA.
- ISO 11133:2014. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- ISO 11290-1:2017 Standard. Microbiology of the food chain. Horizontal method for the detection and enumeration of *Listeria monocytogenes* and for *Listeria* spp.- Part 1: Detection Method
- ISO 11290-2:2017 Standard. Microbiology of the food chain. Horizontal method for the detection and enumeration of *Listeria monocytogenes* and for *Listeria* spp.- Part 2: Enumeration Method
- Jantzen, M.M., J. Navas, M. de Paz, B. Rodriguez, W.P. da Silva & M. Nuñez (2006) Evaluation of ALOA plating medium for its suitability to recover high pressure-injured *Listeria monocytogenes* from ground chicken meat. Letters Appl. Microbiol 43:313-317
- Manafi, M. W. Kneifel & S. Bascomb (1991) Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol Rev. 55:3:335-348
- Ottaviani, F., M. Ottaviani & M. Agosti (1997) Esperienza su un agar selettivo e differenziale per *Listeria monocytogenes*. Industrie Alimentari 36:1-3
- Victor Lachica, R. (1990) Selective plating medium for quantitative recovery of food-borne *Listeria monocytogenes*. Appl. Environ. Microbiol. 56:1:167-169
- Watkins, J. & K.P. Sleath (1981) Isolation and enumeration of *Listeria monocytogenes* from sewage, sewage sludge and river water. J. Appl. Bacteriol. 50:1-9
- UNE-EN ISO 11133 (2014). Microbiología de los alimentos para consumo humano, alimentación animal y agua.- Preparación, producción, conservación y ensayos de rendimiento de los medios de cultivo.

STORAGE

For laboratory use only. Keep tightly closed, away from bright light, in a cool dry place (+4 °C to 30 °C).

