

TECHNICAL DATA SHEET

Article No. 8267

R2A Agar Ph. Eur.

SYNONYMS

-

SPECIFICATION

Solid medium for the enumeration of heterotrophic microorganisms in treated waters, according to Pharmacopoeial Methods.

FORMULA* IN G/L

Proteose peptone	0.500
Casein hydrolysate (Tryptone)	0.500
Yeast extract	0.500
D(+)-Glucose	0.500
Starch	0.500
Sodium pyruvate	0.300
Dipotassium phosphate	0.300
Magnesium sulphate (anhydrous)	0.024
Agar	15.000

Final pH 7.2 ±0.2 at 25 °C

*Adjusted and/or supplemented as required to meet performance criteria.

DIRECTIONS

Suspend 18.1 g of powder in 1 l of distilled water and heat to boiling. Distribute into suitable containers and sterilize by autoclaving at 121 °C for 15 minutes.

DESCRIPTION

R2A Agar was proposed in 1979 by Reasoner and Geldenreich and a few years later accepted by the APHA as an alternative medium for the enumeration of stressed cells in treated potable water.

The use of nutrient rich media like PCA or TSA allows the growth of most microbes, but does not permit the recuperation of stressed or chlorine resistant organisms. Using a medium like R2A with low nutrients in

Th. Geyer GmbH & Co. KG

Dornierstr. 4 – 6
D-71272 Renningen
Tel.: +49 7159 1637-0
Fax: +49 7159 1637-710
renningen@thgeyer.de
www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600)
IBAN DE85600501010002036302
Postbank Stuttgart (Swift/BIC PBNKDEFFXXX)
IBAN DE3260010070000020708
Deutsche Bank (Swift/BIC DEUTDESSXXX)
IBAN DE06600700700125518100

St.-Nr. 70093/40018 / USt-IdNr. DE147510304
Amtsgericht Stuttgart / HRA-Nr. 254140
Persönlich haftende Gesellschafterin:
Geyer Beteiligungsgesellschaft mbH
Amtsgericht Stuttgart / HRB-Nr. 252035
Geschäftsführer: Lutz-Alexander Geyer / Thomas Roth

combination with a lower temperature and longer incubation time it is possible to induce the resuscitation of these damaged cells.

In R2A Agar the source of nitrogen is the peptone and Yeast Extract supplies the vitamins and growth factors. The source of carbon is dextrose and magnesium sulfate and potassium phosphate maintain the osmotic pressure. The starch is a detoxifier and sodium pyruvate increases the recuperation of stressed cells. The agar acts as gelling agent.

TECHNIQUE

The water sample must be processed as quickly as possible. If it is not possible to process within the first 6 hours, the sample must be refrigerated, but not for more than 30 hours.

R2A Agar can be used for pour plates, streak plates or filtration. The pour plate method can affect the recovery capacity of the medium because due to thermal shock when mixing molten agar with the sample. Incubating at 35 °C, an incubation period of 3-5 days is recommended. In most circumstances an incubation temperature of 20-28 °C for 5-7 days is more effective. Plates must be protected against dehydration.

CAUTIONS AND LIMITATIONS

None

QUALITY CONTROL

- Incubation temperature: 30-35 °C/20-28 °C
- Incubation time: 48-72 h/5-7 d
- Inoculum: Practical range 50-100 CFU (productivity), according to Ph. Eur. And ISO 11133:2014. Membrane filter method.

Microorganism	Growth	Remarks
<i>Bacillus subtilis</i> ATCC® 6633	Productivity >0.70	None
<i>Staphylococcus aureus</i> ATCC® 6538	Productivity >0.70	None
<i>Pseudomonas aeruginosa</i> ATCC® 9027	Productivity >0.70	None
<i>Escherichia coli</i> ATCC® 25922	Productivity >0.70	None
<i>Salmonella abony</i> NCTC® 6017	Productivity >0.70	None
<i>Candida albicans</i> ATCC® 10231	Productivity >0.70	None
<i>Aspergillus brasiliensis</i> ATCC® 16404	Productivity >0.70	None

REFERENCES

- ATLAS, R.M. (1995) Handbook of Media for Environmental Microbiology. CRC Press. Boca Raton. Fla. USA.

Th. Geyer GmbH & Co. KG

Dornierstr. 4 – 6
D-71272 Renningen
Tel.: +49 7159 1637-0
Fax: +49 7159 1637-710
renningen@thgeyer.de
www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600)
IBAN DE85600501010002036302
Postbank Stuttgart (Swift/BIC PBNKDEFFXXX)
IBAN DE32600100700000020708
Deutsche Bank (Swift/BIC DEUTDESSXXX)
IBAN DE06600700700125518100

St.-Nr. 70093/40018 / USt-IdNr. DE147510304
Amtsgericht Stuttgart / HRA-Nr. 254140
Persönlich haftende Gesellschafterin:
Geyer Beteiligungsgesellschaft mbH
Amtsgericht Stuttgart / HRB-Nr. 252035
Geschäftsführer: Lutz-Alexander Geyer / Thomas Roth

- CLESCERI, L.S., A.E. GREENBERG and A.D. EATON (1998) Standard Methods for the Examination of Water and Wastewater. 20th ed. APHA Washington D.C. USA.
- EATON, A.D., A.E. GREENBERG and L.S. CLESCERI (1995). Standard Methods for the Examination of Water and Wastewater. 19th ed. APHA Washington D.C. USA.
- EUROPEAN PHARMACOPOEIA. 6th ed. Suppl 6.3 (2009) General Monographs. Water for injections. (pg. 4339) EDQM. Council of Europe. Strasbourg.
- GREENBERG, A.E., R.R. TRUSSELL and L.S. CLESCERI (1985). Standard Methods for the Examination of Water and Wastewater. 16th ed. APHA-AWWA-WPCF. Washington D.C. USA.
- REASONER, D.J. and E.E. GELDREICH (1979) A new Medium for the enumeration and subculture of bacteria from potable water. Abstracts of Annual Meeting. ASM 79th Meeting. Paper #N7.
- ISO 11133:2014. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- Van SOETSBERGER, A.A. and C.H. LEE (1969) Pour plates or streak plates? Appl. Microbiol. 18:1092 -1094.

STORAGE

Keep tightly closed, away from light, in a dry place (4-30 °C).

SHELF LIFE

4 years from date of production.

Th. Geyer GmbH & Co. KG

Dornierstr. 4 – 6
D-71272 Renningen
Tel.: +49 7159 1637-0
Fax: +49 7159 1637-710
renningen@thgeyer.de
www.thgeyer.de

BW-Bank (Swift/BIC SOLADEST600)
IBAN DE85600501010002036302
Postbank Stuttgart (Swift/BIC PBNKDEFFXXX)
IBAN DE32600100700000020708
Deutsche Bank (Swift/BIC DEUTDESSXXX)
IBAN DE06600700700125518100

St.-Nr. 70093/40018 / USt-IdNr. DE147510304
Amtsgericht Stuttgart / HRA-Nr. 254140
Persönlich haftende Gesellschafterin:
Geyer Beteiligungsgesellschaft mbH
Amtsgericht Stuttgart / HRB-Nr. 252035
Geschäftsführer: Lutz-Alexander Geyer / Thomas Roth