

PIPETY BRAND			Transferpette jednokanałowa							Transferpette wielokanałowa						Transferpette S jednokanałowa						Transferpette S wielokanałowa									
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr kategorii	0,1-1	0,5-10	2-20	5-50	10-100	20-200	25-250	100-1000	500-5000	0,5-10	2-20	2,5-25	5-50	10-100	20-200	30-300	0,1-1	0,5-10	2-20	10-100	20-200	30-300	100-1000	0,5-10	5-50	10-100	20-200	30-300	
			BEZ FILTRA	0,5-10	30011 34011	•	•								•								•						•		
30012 34012	•	•											•								•						•				
30013 34013	•	•											•								•						•				
30014 34014	•	•											•								•						•				
30015 34015	•	•											•								•						•				
1-200	30201 30221				•	•	•	•	•				•	•	•	•	•	•			•	•	•				•	•	•	•	•
	30202 30222				•	•	•	•	•				•	•	•	•	•	•			•	•	•				•	•	•	•	•
	30203 30223				•	•	•	•	•				•	•	•	•	•	•			•	•	•				•	•	•	•	•
	30224				•	•	•	•	•	•			•	•	•	•	•	•			•	•	•				•	•	•	•	•
	30225				•	•	•	•	•	•			•	•	•	•	•	•			•	•	•				•	•	•	•	•
1-300	30301						•	•	•						•	•	•				•	•					•	•	•	•	•
	30302						•	•	•						•	•	•				•	•					•	•	•	•	•
	30303						•	•	•						•	•	•				•	•					•	•	•	•	•
100-1000	31001										•															•					
	31002										•															•					
	31003										•															•					
500-5000	35001										•																				
Z FILTREM	0,5-10	30016 34016		•	•								•	•							•						•				
	0,5-20	30026			•								•								•										
	0,5-50	30056				•								•													•				
	1-100	30106				•								•							•						•				
	1-200	30206					•	•						•							•				•			•			
	1-300	30306						•							•												•			•	
	100-1000	31006								•																•					

• kompatybilne

• niekompatybilne

PIPETY CAPP			CappAero Comfort jednokanałowa					CappAero wielokanałowa					ecopipette jednokanałowa					CappTronic jednokanałowa			CappTronic wielokanałowa							
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr kategorii	0,5-10	2-20	5-50	10-100	20-200	100-1000	0,5-10	2-20	5-50	10-100	20-200	30-300	0,2-2	0,5-10	2-20	5-50	10-100	20-200	100-1000	2-20	10-200	100-1000	2-20	10-200	100-1000	
			BEZ FILTRA	0,5-10	30011 34011	●												●	●									
30012 34012	●								●							●	●											
30013 34013	●								●							●	●											
30014 34014	●								●							●	●											
30015 34015	●								●							●	●											
1-200	30201 30221				●	●	●					●	●					●	●	●			●			●		
	30202 30222				●	●	●				●	●	●					●	●	●			●			●		
	30203 30223				●	●	●				●	●	●					●	●	●			●			●		
	30224				●	●	●				●	●	●					●	●	●			●			●		
	30225				●	●	●				●	●	●					●	●	●			●			●		
1-300	30301					●	●					●	●	●					●	●			●			●		
	30302					●	●					●	●	●					●	●			●			●		
	30303					●	●					●	●	●					●	●			●			●		
100-1000	31001							●												●				●			●	
	31002							●												●				●			●	
	31003						●												●				●			●		
Z FILTREM	0,5-10	30016 34016	●												●													
	0,5-20	30026																										
	0,5-50	30056			●												●											
	1-100	30106				●					●							●										
	1-200	30206					●					●							●			●			●			
	1-300	30306											●							●					●			
	100-1000	31006						●												●				●		●		

● kompatybilne

● niekompatybilne

PIPETY EPPENDORF			Reference 2 jednokanałowa						Reference 2 wielokanałowa			Research Plus jednokanałowa						Research Plus wielokanałowa						
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr kategorii	Reference 2 jednokanałowa						Reference 2 wielokanałowa			Research Plus jednokanałowa						Research Plus wielokanałowa						
			0,1-2,5	0,5-10	2-20	10-100	30-300	100-1000	0,5-10	10-100	30-300	0,1-2,5	0,5-10	2-20	5-100	10-100	20-200	30-300	100-1000	500-5000	0,5-10	10-100	30-300	
BEZ FILTRA	0,5-10	30011 34011	•	•					•			•	•								•			
		30012 34012	•	•								•	•									•		
		30013 34013	•	•								•	•									•		
		30014 34014	•	•								•	•									•		
		30015 34015	•	•								•	•									•		
	1-200	30201 30221			•	•	•			•	•			•		•	•	•					•	•
		30202 30222			•	•								•		•	•						•	•
		30203 30223			•	•								•		•	•						•	•
		30224			•	•								•		•	•						•	•
		30225			•	•								•		•	•						•	•
	1-300	30301				•	•			•	•					•	•	•					•	•
		30302				•				•						•	•						•	•
		30303				•										•	•						•	•
	100-1000	31001																						•
		31002																						•
		31003																						•
	500-5000	35001																						•
	Z FILTERM	0,5-10	30016 34016	•	•					•			•									•		
		0,5-20	30026			•							•											
0,5-50		30056																						
1-100		30106				•			•						•							•		
1-200		30206					•			•						•	•					•	•	
1-300		30306																						
100-1000		31006																					•	

• kompatybilne

• niekompatybilne

• kompatybilne, ale wyrzutnik nie  
zrzuca końcówki

PIPETY GILSON			Pipetman Classic jednokanałowa						Pipetman G jednokanałowa						Pipetman L jednokanałowa					Pipetman L wielokanałowa				Pipetman Neo jednokanałowa					Pipetman Neo wielokanałowa			
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr kategorii	0,2-2	1-10	2-20	20-100	50-200	200-1000	0,2-2	1-10	2-20	10-100	20-200	100-1000	0,2-2	1-10	2-20	10-100	20-200	100-1000	0,5-10	2-20	20-200	20-300	0,2-2	1-10	2-20	10-100	20-200	100-1000	2-20	20-200
			BEZ FILTRA	0,5-10	30011 34011	●	●					●	●					●	●								●	●				
30012 34012	●	●								●	●					●	●								●	●						
30013 34013	●	●								●	●					●	●								●	●						
30014 34014	●	●								●	●					●	●								●	●						
30015 34015	●	●								●	●					●	●								●	●						
1-200	30201 30221				●	●	●					●	●	●				●	●	●			●			●	●	●		●	●	
	30202 30222				●	●	●					●	●	●				●	●	●			●			●	●	●		●	●	
	30203 30223				●	●	●					●	●	●				●	●	●			●			●	●	●		●	●	
	30224				●	●	●					●	●	●				●	●	●			●			●	●	●		●	●	
	30225				●	●	●					●	●	●				●	●	●			●			●	●	●		●	●	
1-300	30301					●	●					●	●					●	●				●			●	●		●	●		
	30302					●	●					●	●					●	●				●			●	●		●	●		
	30303					●	●					●	●					●	●				●			●	●		●	●		
100-1000	31001							●												●									●			
	31002							●												●									●			
	31003						●												●									●				
Z FILTREM	0,5-10	30016 34016	●	●					●	●					●	●								●	●							
	0,5-20	30026			●						●						●												●			
	0,5-50	30056																														
	1-100	30106				●						●						●								●						
	1-200	30206					●						●						●							●		●	●			
	1-300	30306																														
	100-1000	31006					●												●									●				

● kompatybilne

● niekompatybilne

PIPETY SOCOREX			Accura Manual jednokanałowe							Accura Manual wielokanałowe				Accura Manual XS jednokanałowe										
KOŃCÓWKI HTL	Zakres pojemności końcówek [μl]	Nr kategorii	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000	0,5-10	5-50	20-200	40-350	0,1-2	0,5-10	1-10	2-20	5-50	10-100	20-200	100-1000		
BEZ FILTRA	0,5-10	30011 31011	●	●	●						●				●	●	●							
		30012 34012	●	●	●							●				●	●	●						
		30013 34013	●	●	●							●				●	●	●						
		30014 34014	●	●	●							●				●	●	●						
		30015 34015	●	●	●							●				●	●	●						
	1-200	30201 30221				●		●	●				●	●	●				●		●	●		
		30202 30222				●		●	●				●	●	●				●		●	●		
		30203 30223				●		●	●				●	●	●				●		●	●		
		30224				●		●	●				●	●	●				●		●	●		
		30225				●		●	●				●	●	●				●		●	●		
	1-300	30301						●	●				●	●	●					●	●	●		
		30302						●	●				●	●	●					●	●	●		
		30303						●	●				●	●	●					●	●	●		
	100-1000	31001									●													●
		31002									●													●
31003										●													●	
Z FILTREM	0,5-10	30016 34016	●	●	●						●				●	●	●							
	0,5-20	30026			●	●												●						
	0,5-50	30056					●					●							●					
	1-100	30106						●												●				
	1-200	30206							●				●	●								●		
	1-300	30306																						
	100-1000	31006								●													●	

● kompatybilne

● niekompatybilne

PIPETY THERMO			Finnpipette F1 jednokanałowa							Finnpipette F1 wielokanałowa				Finnpipette F2 jednokanałowa							Finnpipette F2 wielokanałowa				Finnpipette F3 jednokanałowa						Finnpipette F3 wielokanałowa							
KOŃCÓWKI HTL	Zakres pojemności końcówek [µl]	Nr kategorii	0,2-2	0,5-5	1-10	2-20	5-50	10-100	20-200	30-300	100-1000	1-10	5-50	10-100	30-300	0,2-2	0,5-5	1-10	2-20	10-100	20-200	30-300	100-1000	500-5000	1-10	5-50	10-100	30-300	0,2-2	1-10	5-50	2-20	10-100	20-200	100-1000	5-50	30-300	
			BEZ FILTRA	0,5-10	30011 34011	•	•	•							•			•	•	•							•				•	•						
30012 34012	•												•			•									•				•									
30013 34013	•												•			•									•				•									
30014 34014	•												•			•									•				•									
30015 34015	•												•			•									•				•									
1-200	30201 30221								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	30202 30222								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	30203 30223								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	30224								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
1-300	30225								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	30301								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	30302								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
100-1000	30303								•	•	•			•	•	•					•	•	•			•	•	•		•					•	•		
	31001											•																									•	
	31002											•																									•	
500-5000	31003										•																										•	
	35001																																					•
Z FILTREM	0,5-10	30016 34016	•	•	•							•			•	•	•							•				•	•									
	0,5-20	30026											•													•												
	0,5-50	30056											•													•										•		
	1-100	30106						•											•										•									
	1-200	30206							•	•																									•			
	1-300	30306								•				•														•									•	
	100-1000	31006									•																										•	

• kompatybilne

• niekompatybilne