

## Dual-Action Shaker KL 2



### Description

Depending on the task the shaking motion of the Dual Action Shaker KL 2 can be switched from orbital to reciprocating.

### Basic Equipment

Basic device incl. rack system with 3 clamping strips and rubber mat

## Loading Capacity

Erlenmeyer flasks 250 ml 9 pieces

Separating funnels 250 ml 1 pieces

You will find further loading possibilities [here](#)

## Details

**Motion:** can be switched from orbital to reciprocating (horizontal)

**Shaking platform:** 300 x 300 mm

**Max. load:** 2 kg

**Shaking speed:** 30 - 420 rpm, steplessly adjustable

**Stroke:** 8 mm

**Timer:** 0 - 120 min / continuous

**Order No.:** 6115 000

**Accessory:** Clamping strip h 0050 118

Clamping strip v 0050 477

Cramp strip 0050 206

Spring strip 0050 207

Test tube rack with 44 drillings à 14 mm Ø

Test tube rack with 44 drillings à 16 mm Ø

Test tube rack with 44 drillings à 18 mm Ø

Test tube rack with 14 drillings à 30 mm Ø

## Technical specifications

**Electrical supply:** 230 V or 115 V, 50 / 60 Hz please indicate in case of order

**Enclosure protection:** IP 21

**Heat emission:** approx. 4 W

**Ambient temperature:** 5°C bis 50°C

**Relative humidity:** ~ 85 %

**Dimensions (w x d x h):** 355 x 455 x 195 mm

**Weight:** 13,5 kg