

TECHNICAL DATA SHEET

Article No. 9360

LB Broth acc. to Miller, ready-to-use culture medium

SPECIFICATION

Prepared medium. Liquid medium for general purposes and molecular genetics studies of *Escherichia coli*.

Colour: Straw-coloured yellow
pH: 7.0 ± 0.2 at 25°C

COMPOSITION IN G/ L

Casein peptone	10.00
Sodium chloride	10.00
Yeast extract	5.00

PACKAGE DETAILS

9360-10x100ML

Volume	100 ± 3 ml
Bottle size	125 ml
Packaging unit	10 bottles

1 box with 10 x 100 ml in 125-ml-bottles. Injectable cap; Plastic screw inner cap.

DESCRIPTION/TECHNIQUE

Description

LB Medium was originally formulated by Luria and Burrous, Lennox added sodium chloride to improve the osmolarity of the media. The formulation of this broth is according to the Lennox recipe, modified by Miller, who has increased the sodium chloride concentration.

Technique

Collect, dilute and prepare isolates and cultures as required according to specifications, directives, official standard regulations and/or expected results. Dispense liquid medium in appropriate containers if the original container is of large volume.

Inoculate aseptically the bottles with the prepared sample or its dilution. Incubate the bottles tightly closed aerobically at 35-37°C for 24-48 h. Incubation times, temperature and sample volumes may vary depending on the sample, on the specifications.



This medium can be used to inoculate any confirmatory, secondary medium by streaking methodology or by spiral method. After incubation, enumerate all the colonies that have appeared onto the surface of the secondary agar. Each laboratory must evaluate the results according to their specifications.

MICROBIOLOGICAL CONTROL

Prepare tubes - Inoculate: Practical range 100 ± 20 CFU. min. 50 CFU (productivity).

Microbiological control according to the current version of the ISO 11133:2014/A1:2018.

Analytical methodology according to ISO 11133:2014/A1:2018; A2:2020.

Aerobiosis. Incubation at 35 ± 2 °C, reading after 18-24 hours.

Microorganism	Growth
<i>Escherichia coli</i> ATCC® 25922, WDCM 00013	Good
<i>Escherichia coli</i> ATCC® 8739, WDCM 00012	Good
<i>Escherichia coli</i> ATCC® 11775, WDCM 00090	Good
<i>Escherichia coli</i> ATCC® 35218	Good

Sterility control:

Incubation 48 hours at 30-35 °C and 48 hours at 20-25 °C: NO GROWTH.

Check at 7 days after incubation in same conditions.

REFERENCES

- ATLAS, R.M. & L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- AUSUBEL, F.M., R. BRENT, R.E. KINGSTON, D.D. MOORE, J.G. SEIDMAN, J.A. SMITH & K. STRUHL (1994) Current protocols in molecular biology. Greene Pub. Assoc. Inc. Brooklyn. NY.
- GHERNA, R., P. PIENTA & R. COTE (Eds.) 1992. ATCC Catalogue of Bacteria and Bacteriophages. Media #1065, #1226, #1226,
- #1235, #1236, #1315, #1364. American Type Culture Collection. Rockville. MD. USA.
- ISO 11133:2014/ Adm 1:2018. Microbiology of food, animal feed and water. Preparation, production, storage and performance testing of culture media.
- LURIA, S.E. & J.W. BURROUS (1955) Hybridization between *Escherichia coli* and *Shigella*. J. Bacteriol. 74:461-476.
- LENNOX, E.S. (1955) Transduction of linked genetic character of the host bacteriophage P1. Virology 1:190-206.
- MILLER, J.H. (1972) Experiments in Molecular Genetics. Cold Spring Harbor Laboratory. Cold Spring Harbor. NY.
- SAMBROOK, J., E.F. FITSCH & T. MANIATIS (1989) Molecular cloning: A laboratory manual. 2nd ed. Cold Spring Harbor Laboratory. Cold Spring Harbor. NY.
- MILLER, J.H. (1992) A short course in bacterial genetics: A laboratory manual and handbook for *Escherichia coli* and related bacteria, p. 194-195. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, N.Y.

STORAGE

8 - 25 °C

SHELF LIFE

16 months unopened from date of manufacture

updated: 05.06.2023

